

Fero wooden a-sync

16425 The Fero relax chair with castors and lowerable armrests is equipped with a relax mechanism that has been especially developed for the care sector. The a-sync mechanism allows the backrest and legrest to be adjusted separately.

Support

- 4 beech laminated legs, connected to each other by means of solid beech support beams.
- Glued mortise and tenon joint.
- The sides are firmly connected to each other by means of a solid beech cross beam and a metal frame.
- Comfortable rounded armrests of $\pm 7,5$ cm.
- All edges and corners have been rounded.
- 4 Integral S design castors \varnothing 125 mm on precision ball bearings, 2 castors at the back have central brakes with a directional lock.
- Extending, synthetic footrest with a step safety and with an antiskid structure.

Seat

- Removable seat, made of a melaminated MDF board of 12 mm.
- Upholstered with fire-retardant foam in 2 different density zones to guarantee an excellent seat comfort.
- The rounded seat in front stimulates a fluent blood circulation.
- Comfortable seat height of 47,5 cm.
- The upholstery is finished with imitation leather and decorative stitching.

Backrest

- Made of an ergonomically preformed plywood shell of 12 mm.
- Upholstered with fire-retardant foam in 3 different density zones to guarantee an excellent seat comfort.
- Equipped with an adjustable and removable head/neck cushion.
- Push bar at the back.

Legrest

- Made of a plywood board of 10 mm, upholstered with fire-retardant foam and finished with imitation leather.

Lowerable armrests

- The ergonomic front design of the beech laminated armrests makes it easy to stand up or sit down.
- Standard equipped with dinner tray attachments.
- Two solid chromium-plated axes to adjust the armrests in height in 4 different positions.
- Locked downwards by means of a spring lever.
- In their lowest position, the armrests are on the same level as the seat.

A-sync mechanism

- Allows the backrest and the legrest to be adjusted separately.
- The backrest can be adjusted by means of a lever on the right side of the chair.
- A gas spring minimizes the necessary activation strength of the mechanism and ensures a continuous adjustment of the backrest.
- The legrest and seat can be easily adjusted by just pushing yourself off a little, without needing a lever.
- Optimal leg support, thanks to the telescopic legrest.
- When the chair is put in the seat position, the legrest is turned inwards so that you can stand up without any problem.

Colours

We would like to refer to our extensive colour gamut.
Sample card + samples available on request.

Accessories


We would like to refer to our
accessorielist.



Dimensions

16425 | Fero wooden relax chair:

- L 68 x D 65-155 x H 120,5 cm
- Seat height: 47,5 cm
- Seat depth: 51 cm
- Seat width: 53,5 cm
- Height armrests: 65,5 cm
- Safe working load: 135 kg



Materials and finishing

- Wood: beech (solid, laminated), Plywood, MDF
- Metal: epoxy coating (light grey).
- Fire-retardant foam.
- Upholstery: see collection.
- Resistant to the common cleaning products.

For further information about used **materials, constructions and maintenance**: consult our materialfiles.

Options

- Elice dinner tray.
- Abduction block.
- Lock mechanism for the lever.
- Seat cushion with visco-elastic foam and a bi-elastic cover.
- Anatomical reductor.
- Accessory brackets.
- Restraint bracket.
- Lateral headrest.
- Case for IV rod.
- Backrest with removable headrest.
- Dinner tray support on backrest.
- Double design brake pedal.
- Twin castors \varnothing 125 mm.

